

AED



AUTOMATIC EXTERNAL DEFIBRILLATOR

Accurate and Reliable ECG Analysis
Easy to Use
Motion Detector
Pace Maker Detection

Self Test Function
AHA/ERC Guidelines for CPR and Defibrillation
5 Year Warranty



Turkish First Aid Association supports "Life-Point"

MedicinskáTechnika.CZ
Internetový obchod se zdravotnickou technikou
servis a prodej zdravotnické techniky

+420 603 440 557 | info@medicinskathechnika.cz | www.medicinskathechnika.cz



**AUTOMATIC
EXTERNAL
DEFIBRILLATOR**

Life-Point Pro AED Technical Specifications

| | |
|------------------------------|--|
| Operating Mode | Automatic |
| Charging Time at max Joule | typically < 10 sec |
| ECG Analysis Period | typically < 10 sec |
| Wave Form | Truncated Biphasic Wave Form |
| Energy Steps | Adult : 150J - 150J - 200J / Child : 50J - 50J - 50J |
| Pace Maker Detection | Yes |
| Motion Detector | Yes |
| Impedance Measurement | 25 – 175 ohm |
| Maximum Discharging | typically 100 shocks (at 20 °C with full battery) |
| Derivation | D II |
| ECG Bandwidth | 2 – 25 Hz |
| Battery Type | LiMnO ₂ (12V 4500 mA) |
| Self Test Schedule | Daily - Weekly - Monthly - Each Start |
| Self Test Scope | Battery - Internal Electronics - Software |
| Dimensions | 9,5cm x 27cm x 24cm |
| Weight (battery included) | 1850 gr |
| Class (93/42/EEC) | II B |
| Operating Conditions | 0°C - 50°C ve % 10 – % 95 RH |
| Storage Conditions | Between -20°C and 60°C |
| Data Documentation | ECG Signal/ Event / Voice Record |
| Data Transfer | Metsis AED Software (mini USB) |
| Minimum Computer Requirement | Windows XP – Windows Vista – Windows 7 Pentium or higher VGA Monitor or higher CD-ROM Driver USB Port 2 Gb Free Disk Space |

www.metsismedikal.com



Merkez | Head Office
İvedik Organize Sanayi Bölgesi
Soğutmaçlar Sit. 1444. Sokak No: 11
Ostim - Ankara / Turkey
Tel: +90.312 395 77 70
Faks: +90.312 395 77 72

Ar-Ge | R&D
ODTÜ Ostim Teknokent İdari Binası
Uzay Çığı Caddesi 78. Sokak No: 6
2.B Kat D: 1-2 Ostim-Ankara/Turkey
Tel: +90.312 386 21 90
Faks: +90.312 386 21 80

info@metsismedikal.com